



## CircuFlow™ 5150 Intermittent Pneumatic Compression (IPC) Pump

The CircuFlow™ 5150 Intermittent Pneumatic Compression (IPC) Pump is used to treat chronic lymphatic and venous disorders of the lower extremities and chronic lymphatic disorders of the upper extremities.



### Key Features:

- FDA and CE Approved
- Use indicator timer which shows patient use and compliance
- Quiet, wave-like gradient sequential inflation/deflation cycle provides patient with comfortable massage-like therapy.
- Easy-to-read LCD screen that allows the user to see real-time pressure in each individual chamber.
- Convenient digital controls, built-in timer feature for easy therapy time control.
- Automatic shut-off and emergency stop feature for maximum user safety assurance.
- Vacuum mechanism evacuates air from the garments upon shut-off of the system allowing for greater ease in the applying and taking off the garments.
- Quick-release hose connectors allow for easy removal of the garments for cleaning as needed.

### Specifications

- Dimension:** 10" (W) x 9.5" (D) x 4.5" (H)
- Weight:** 5.5 lbs
- Inflation Cycle:** 60 seconds
- Rest Time:** 12 seconds
- Cycle Time:** 15 seconds / chamber
- End Therapy Vacuum Time:** 180 seconds
- Treatment Time:** 0.5 – 3 Hours
- Pressure Range:** 20 – 80 mmHg
- Sleeve Prefill:** To 10 mmHg
- Gradient:** 7%; example 3 mmHg at 50 mmHg Distal Pressure
- Electrical:** 120VAC, 60Hz, 45 VA maximum
- Fuse Rated:** 250VAC, 1.0 AMP, SLO-BLO
- Protection Against Electric Shock:** Class II
- Operation Mode:** Intermittent Gradient Sequential

LB034.0002 Rev C 20140312



## CircuFlow™ 4-Chamber Garment Measurement Chart

Lower Extremity		Model	Length (in.)		Circumference (in.)		
Size	Type	4 chamber	Inseam	Ankle	Knee	Thigh/Groin**	
Small	Full Leg	D-300S	23	21	25	30	
Small Wide	Full Leg	D-300SW	26	24	29	35	
Small Extra Wide	Full Leg	D-300SXW	24	28	35	40	
Medium	Full Leg	D-300M	29	22	27	31	
Medium Wide	Full Leg	D-300MW	29	25	31	35	
Medium Extra Wide	Full Leg	D-300MXW	28	29	36	40	
Large	Full Leg	D-300L	34	21	27	31	
Large Wide	Full Leg	D-300LW	36	22	31	40	
Half	Half Leg	D-301H	16.5	22	28	-	
Half Wide	Half Leg	D-301HW	19	24	30	-	

  

Upper Extremity		Model	Length* (in.)		Circumference (in.)		
Size	Type	4 chamber	Axilla to Fingertips	Axilla	Elbow	Wrist	
Medium	Arm	D-302M	26.5	22	19	13	
Large	Arm	D-302L	32	23	20	14	

Lower Extremity		Model	Length (cm.)		Circumference (cm.)		
Size	Type	4 chamber	Inseam	Ankle	Knee	Thigh/Groin**	
Small	Full Leg	D-300S	58.4	53.3	63.5	76.2	
Small Wide	Full Leg	D-300SW	66	61	73.7	88.9	
Small Extra Wide	Full Leg	D-300SXW	61	71.1	88.9	101.6	
Medium	Full Leg	D-300M	73.7	55.9	68.6	78.8	
Medium Wide	Full Leg	D-300MW	73.7	63.5	78.7	88.9	
Medium Extra Wide	Full Leg	D-300MXW	71.1	73.7	91.4	101.6	
Large	Full Leg	D-300L	86.4	53.3	68.6	78.7	
Large Wide	Full Leg	D-300LW	91.4	55.9	78.7	101.6	
Half	Half Leg	D-301H	41.9	55.9	71.1	-	
Half Wide	Half Leg	D-301HW	48.3	61	76.2	-	

  

Upper Extremity		Model	Length* (cm.)		Circumference (cm.)		
Size	Type	4 chamber	Axilla to Fingertips	Axilla	Elbow	Wrist	
Medium	Arm	D-302M	67.3	55.9	48.3	33	
Large	Arm	D-302L	81.3	58.4	50.8	35.6	

\* Length corresponds to the entire sleeve length from axilla to fingertips \*\* The circumference of the upper rim of the sleeve

LB18.0001 Rev C 20140312