

Lympha Press Optimal®



MODEL	Lympha Press Optimal® I201-AP
CATALOG NUMBER	LI0003C6
CYCLE	<ul style="list-style-type: none"> • Pretherapy™ • Sequential • Peristaltic
PRESSURE RANGE	20-90 mmHg
WEIGHT	6.5 kg / 14.3 lbs
SIZE	15x32x38.5 cm / 6x12.5x15.3 in
ELECTRICAL	115V / 60 Hz
ACCESSORIES	Hose Bundle L25S1000



GARMENT COMPATIBILITY

Lympha Press Optimal® is used with our 12 to 24 overlapping chamber appliances, including leg and arm garments, Lympha Pants™, Comfysleeve™ I-75, and Lympha Jacket II™. For use with LymphaPod®, select Optimal Model I201-AP-LPOD.

Lympha Press Optimal® is HCPCS coded E0652.

PLEASE SEE PAGE 32 FOR INDICATIONS/CONTRAINDICATIONS.

Advanced pneumatic compression with Pretherapy™.

Lympha Press Optimal® combines the gentle touch of manual lymph drainage with Lympha Press® calibrated gradient compression in a home care system that is easy to use.

Optimal® offers three different treatment modes, including: Pretherapy™ (based on the principles of manual lymph drainage), peristaltic, and sequential. Pressure can be adjusted independently in each chamber or in four zones of three chambers each to reduce pressure over sensitive areas. Therapists appreciate that pressures and modes are locked in for patient compliance. The clear digital display makes it easy to follow the progress of each treatment session. And the Lympha Press Optimal® has another unique feature – pressures and modes can be set using our software and a USB connection.

FEATURES

- Choice of Pretherapy™, peristaltic, and sequential modes.
- 12-zone gradient calibrated “at the pump.” Pressure can be adjusted in each compression chamber individually, or in four groups of three chambers each.
- Garments for torso clearance. Use with Lympha Pants™, Lympha Jacket™, or ComfySleeve™ I-75 for torso treatment on the principles of manual lymph drainage.
- Bilateral or unilateral treatment.
- Patient lock – pressure and modes are set and locked, according to prescription.